



Clinical Swing-Out Centrifuge (iFuge CS03)

Clinical Swing-Out Centrifuge




- Unique 3 tube rotor design (120 deg. balanced rotor) - allows use of 50 ml swing in a smaller compact unit
- Maintenance free, long life BLDC Motor
- Speed Selection: 500 - 3200 RPM (in steps of 10RPM) – Delivers upto 1470 x g
- Timer setting for 1 to 999 mins
- Smart airflow design ensures chamber temperature rise is under control
- Imbalance detection feature with Auto-cut off safety

Clinical Swing-Out Centrifuge

- Digital display with speed & time setting
- Favorite mode: Allows user to set upto 9 User Defined Programs
- Lid lock safety feature – lid opens automatically on run completion
- Low loading height for ease of operation.
- Swing out 3 x 50 ml tube rotor with adaptors for 9,6,4,3 & 2ml tubes
- One touch RPM / RCF conversion & short spin operation option on keypad



Advantages of Clinical Swing-Out Centrifuge



- A compact unit best suited for applications that needs clear line separations for various applications for tubes upto 50 ml capacity.

Advantages of 'Horizontal Separation' are:

i) Clear Straight-line separation as shown in the adjacent picture

ii) Tighter separation line eliminates Re-Spins

- Few Applications: Blood (STAT), Urine, Semen, Fecal, PRP, PPP, Stem-Cell etc.



neuation

TECHNOLOGIES

www.neuation.com